PRESENTATION:

1 Syringe with 2.5 mL and 2 tips.

DESCRIPTION:

10% acid gel is an hydrofluoric acid-based gel, with a concentration of 10%, and is the most efficient and used way to attack porcelain surfaces, allowing for a greater and better retention of resinous material of feldspathic porcelain, ceramic glass, and orthodontic braces. Its appropriate viscosity and the presence of dye allow for safe application, without the danger of flow and better visualization of the region to be conditioned.

INDICATION:

Surface conditioning of dental porcelains for adhesion in restorative processes, with feldspathic porcelains and ceramized glasses.

COMPOSITION:

10% hydrofluoric acid, thickening agent, dye, and purified water.

MODE OF USE:**For feldspathic porcelain:**

Use 10% Acid Gel for 1-2 minutes.

For manufactured ceramic glass:

Use 10% Acid Gel for 1 minute.

For pressed/injected ceramic glass:

Use 10% Acid Gel for 20-30 seconds.

For all types of materials it is necessary to use an abundant amount of water to remove the product, since contact between hydrofluoric acid and enamel can cause bone decalcification.

In orthodontic procedures such as installation of braces on porcelain pieces, it is necessary to remove the product with cotton or a surgical.

WARNING AND PRECAUTIONS:

- Restricted use, for professionals who specialize in dentistry, in accordance with local legislation and with specific training for its handling;
- Rigorously follow the procedures provided in the method for use;
- The use of this product for the conditioning of enamel and dentin is strictly prohibited;
- Use gloves and safety glasses while handling product;
- Protect the patient by using safety glasses;
- Do not use the product after it has passed its expiration date;
- Keep out of reach of children.

CONTRAINDICATIONS AND SIDE EFFECTS:

Do not use on patients with sensitivity to the component parts of this product.

STORAGE AND TRANSPORT:

Store in a dry place, away from heat and humidity, in the temperature range between 5 °C and 25 °C. Humidity: below 75 UR.

DISPOSAL:

The product should be disposed of in accordance with active legislation in your country.

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**ÁCIDO GEL
FLUORHÍDRICO 10%
10% HYDROFLUORIC
ACID GEL**

MADE IN BRAZIL

maquira®
Dental Group**Fabricante | Manufacturer:**

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